

Leisure COLLECTION




shown in
Autumn


Autumn




Cornflower



Spring



Rosewood



Chestnut



Burgundy



This traditional style offers a plush, button back with a sleek armrest design that will subtly match many different decorative styles.

UC 214

UltraComfort[™]
AMERICA
Positioning you for life.



FOR A DIGITAL VERSION
OF THIS INFORMATION,
SCAN HERE.



Leisure COLLECTION

UC 214

Model Number	UC214
Size	Medium
Weight Capacity	325 LB
Recline Positions	2
Standard Fabrics	6
Chaise Seat	No
Back Type	Button
Head Rest Cover	Yes
Arm Covers	Yes
Battery Back Up	Yes
Overall Width	31"
Overall Height	41"
Width Between Arms	20"
Floor to Top of Seat	18"
Seat Depth	19"
Seat to Top of Back	27"
Seat to Top of Arm	7"
Seat Angle Raised	16°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	35"
Back Angle Closed	110°
Back Angle Reclined	135°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	160°
Distance Required from Wall Reclined	11"
Extension Reclined from Back of Seat to End of Foot Rest	40"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical parts	1 YR
Pro Rated Years 4 Thru 7	No
Warranty on Mechanical Labor	1 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam	1 YR
Removable Back Option	No
Customizable *	No
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Yes
4-Motor Heat & Massage Available	No
2-Motor Heat & Massage Available	No
Foot Rest Extension Option Available	No
Smart Tek System	Yes

* See Customize Chart for dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



User-friendly feature of every UltraComfort lift and recline chair:

Easy to use battery backup system. Each UltraComfort lift and recliner has a battery backup system in the electrical transformer which operates on two 9 volt batteries. In the event of a power failure, the battery system will raise the chair from the reclined to the seated position. The batteries should be replaced after one use.



UltraComfort™
 AMERICA

Positioning you for life.

UltraComfort
 Phone: 866-793-5100
 Fax: 866-255-9962
 web: www.ultracomfort.net



Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

Listing # E112417

FNAE 80-214A



SHT-UC214.052013